

50 Years of Innovation



ThinProfile™ Platinum



ThinProfile™ Gold



Contour™ Platinum



Contour™ Gold

MedDev Sterile Eyelid Implants

ThinProfile™ Platinum

Standard Sizes — Thickness 0.65 mm / Width 6.0 mm (6.6 mm)*

LL 6006	0.6 grams
LL 6008	0.8 grams
LL 6010	1.0 grams
LL 6012	1.2 grams
LL 6014*	1.4 grams
LL 6016*	1.6 grams
LL 6018*	1.8 grams

ThinProfile™ Gold

Standard Sizes — Thickness 0.65 mm / Width 6.0 mm (6.6 mm)*

LL 4006	0.6 grams
LL 4008	0.8 grams
LL 4010	1.0 grams
LL 4012	1.2 grams
LL 4014*	1.4 grams
LL 4016*	1.6 grams

Contour™ Platinum

Standard Sizes — Thickness 1.0 mm / Width 5.0 mm

LL 5006	0.6 grams
LL 5008	0.8 grams
LL 5010	1.0 grams
LL 5012	1.2 grams
LL 5014	1.4 grams
LL 5016	1.6 grams
LL 5018	1.8 grams

Contour™ Gold

Standard Sizes — Thickness 1.0 mm / Width 5.0 mm

LL 3006	0.6 grams
LL 3008	0.8 grams
LL 3010	1.0 grams
LL 3012	1.2 grams
LL 3014	1.4 grams
LL 3016	1.6 grams
LL 3018	1.8 grams
LL 3020	2.0 grams
LL 3022	2.2 grams
LL 3024	2.4 grams
LL 3026	2.6 grams
LL 3028	2.8 grams

Product Identification: MedDev eyelid weight implants are marked on the concave side with the weight in grams.

Product Packaging: MedDev eyelid weight implants are provided STERILE in individually sealed, double-peel pouches.

MRI Safety: All MedDev Gold Contour™ and ThinProfile™ eyelid weight implants and Platinum Contour™ and ThinProfile™ eyelid weight implants are MR CONDITIONAL. Please refer to 'MedDev MRI Safety Information for MedDev Eyelid Weight Products' (11-10834) available online for additional information.

Special Orders: Special size Contour™ Platinum eyelid weight implants (2.0 grams to 2.8 grams) are available upon request.

ORDERING INFORMATION:

50 Years of Innovation

MedDev External Eyelid Weight Products

Blinkeze® External Lid Weights

Standard Sizes — Thickness 1.0 mm / Width 5.0 mm

TEW 6006	0.6 grams
TEW 6008	0.8 grams
TEW 6010	1.0 grams
TEW 6012	1.2 grams
TEW 6014	1.4 grams
TEW 6016	1.6 grams
TEW 6018	1.8 grams

SKIN TONE SELECTION

Specify the appropriate skin tone suffix when ordering a Blinkeze External Lid Weight: -0 (light), -1 (Pink/Neutral), -2 (Tan), -3 (Medium Brown), -4 (Dark Brown)

LL 4000 Skin Tone Selector Card

Accessories

Tantalum Eyelid Weight Sizing Sets

TEW 7000	Standard Sizing Set 7 weights from 0.6 grams to 1.8 grams
TEW 7004	Special Sizing Set 5 weights from 2.0 grams to 2.8 grams

EyeClose® Double-Sided Adhesive Strips

LL 0505 Box of 100 Strips

Resources

SM 1010	MR Safety Information
SM 1015	Instructional Media
SM 1020	Blinkeze Skin Tone Selector Card/ Physician Order Forms
SM 1025	Clinical References
SM 1030	Billing Guides
SM 1035	Patient Brochures (12 pk.)

Product Identification: MedDev external eyelid weight products are marked on the concave side with the weight in grams.

Product Packaging: MedDev external eyelid weight products are provided NON STERILE. Blinkeze External Lid Weights are packaged in a treatment kit with 100 adhesive strips and an application mirror. Tantalum Sizing Eyelid Weights are packaged in a storage case with 100 adhesive strips and a Skin Tone Selector Card.

MRI Safety: All MedDev external eyelid weight products are MR UNSAFE. Please refer to 'MedDev MRI Safety Information for MedDev Eyelid Weight Products' (11-10834) available online for additional information.

Special Orders: Special size external eyelid weights (2.0 grams to 2.8 grams) are available upon request.

ORDERING INFORMATION:

MedDev
CORPORATION

730 North Pastoria Avenue
Sunnyvale, CA 94085-3522

Telephone:
800.543.2789
408.730.9702

Facsimile:
408.730.9732

Email:
info@meddev-corp.com

Visit us online at www.meddev-corp.com